

Installing Python-Macros in FIJI

Stefan Fischer, Tübingen Structural Microscopy Core Facility (TSM)


August 2024, Version 1.0

The content of this document is shared under Creative Commons ShareAlike 4.0 International (CC BY-SA 4.0)

Although the macros can just be dragged and dropped onto the FIJI Toolbar (as for opening images) to open them, it is much more convenient to install them directly in FIJI so that they are available via the Plugins Menu.

To do so, generate a new folder on your desktop with the corresponding name e.g. “TSM-Macros” and drop all macros of interest into that folder.

Now copy this folder to your FIJI installation under Fiji.app/scripts/Plugins/ and restart Fiji.



The macros can now be directly started via “Plugins” and then from the dropdown list under “TSM-Macros”.